

Cuff for Blood Pressure Monitor



Patented Cone Cuff



Position Guideline



Latex Free



Easy to Apply




Clear Position Guideline



Cone Cuff Patent Nos :
ZL 2003 2 0100815.4 Nr.
20 2004 014 458.1
M244879 Other patents
pending

Storage and Transportation Condition
-10°C~60°C(14°F~140°F) 10%
~90% RH; 700~1060 hPa

			
SMALL	MEDIUM	LARGE	UNIVERSAL
18 - 26 cm	24 - 36 cm	34 - 46 cm	24 - 40 cm
7.0"-10.3"	9.5"-14.2"	13.4"-18.1"	9.4"-15.7"

Distributed by:
Silverline Medical Pty Ltd
8A, 380 Pennant Hills Road, Pennant Hills, NSW 2120
1800 959 777
www.rossmax.com.au