

Product Catalogue Australia & New Zealand 2024

Enteral Nutrition Feeding Tubes, Accessories & Gravity Systems





Freka® GastroTube

Replacement tube for a fully formed and healed stoma



gastric

Gastric Replacement System FR 15 for use within a fully formed and healed stoma. Initial Placement System (PEG) would have been placed at least four weeks previously. Balloon activating lumen connectors will not be changed to ENFit in any way, as they are not within the scope of ISO 80369-3.

- Replacement System of silicone rubber
- X-ray contrast strip in shaft and additional X-ray contrast above balloon (Fig. 1)
- Lateral access with valve for filling balloon (5 - 7.5 ml sterile water or isotonic saline solution)
- Flat, rounded, terminal opening
- Plastic mandrin with guide for Seldinger wire
- Separate ENFit insert for easy connection to the feeding funnel after placement of the feeding tube (Fig. 2)

Freka GastroTube FR 15

FR
15

- For long-term intragastric feeding
- Suitable for all gastric stoma FR 15
- Reliable placement with Seldinger Technique
- Long placement duration with high balloon stability

Contents per single set

- Insert ENFit to convert the Funnel connector to ENFit connector
- Freka GastroTube for FR 15
- Plastic mandrin
- 10 ml Luer syringe

Latex-free
DEHP-free*

*not intentionally made with DEHP

	Sales unit	Art.-no.
Freka GastroTube FR 15	5	7755648